

# MAESTRO<sup>®</sup> PEDICURE SPA

## Specifications

**Dimensions:** 70.4" L (179 cm) x 36.2" W (92 cm) x 55" H (140 cm)

**From Floor to Seat:** 25" (64 cm)

**Shipping Dimensions (Carton Size):** 84" L (214 cm) x 40" W (102 cm) x 61" H (155 cm) The Maestro ships fully assembled.

**Net Weight:** 310 lbs (141 kg)

**Shipping Weight:** 375 lbs (190 kg)

**Client Weight Capacity:** 300 lbs (136 kg)

**Basin Capacity:** 5.9 US Gallons (22.3 liters)

**Power Requirements:** 120V AC, 20 Amp, 60HZ. A **separate and dedicated** 20 Amp GFCI protected circuit **MUST** be used for each foot spa.

**Construction:** Welded steel frame with veneered wood panels. ABS cover panels.

**Basin:** Solid surface, pre-plumbed.

**Jet System:** Two (2) self contained NoPipe™ Jets; removable screen and impeller, no recirculating pipes.

**IMPORTANT: The Maestro Pedicure Spa MUST be installed by a licensed electrician and a licensed plumber in accordance with all national, state and local codes/regulations.**

### Plumbing:

**Incoming:** Hot/Cold supply lines required (1/2" ID), solenoid controlled, built-in fill sensor. Connections located on back of spa or as specified otherwise.

**Sprayer/Filler:** Inline back-flow prevention is provided.

**Discharge:** Gravity drain or optional factory installed drain pump. If factory installed drain pump is utilized, built-in back-flow prevention is provided. Two drain pump sizes available: Standard = 8' Head, High Capacity (optional) = 16' Head

**NOTE:** If through floor installation is desired, refer to Plumbing Rough-in Drawing (available at [ContinuumFootspas.com](http://ContinuumFootspas.com)). See Owners Manual for complete installation instructions.


**Massage:** Rolling and kneading massage, full-back, upper -back or lower-back

**Controls:** Technician's keypad: Seat Forward/Backward, Tilt and Exit. Electronically controlled, low voltage.

**Armchair:** Left Side: Seat Forward/Backward, Tilt and Exit. Electronically controlled, low voltage.

Right Side: Handheld Massage Control Wand

**Chair Release Lever:** Manual operation, locks and unlocks chair to allow swivel.

**Certification:** Underwriters Laboratories Listed  US Compliance Listings to the following Standards: UL 1951, CAN/CSA C22.2 No. 68-92, ANSI/ASME A112.18.3M, A112.19.7M, A112.18.1, CSA B125.1.2005  
Compliant with IAPMO Plumbing Standards. Fixtures and hoses are UPC compliant.

**Limited Warranty\*:** Basin- Sixty (60) months from date of shipment from Continuum Footspas facility.

All other parts and accessories - Twelve (12) months from date of shipment from Continuum Footspas facility.

Labor: Twelve (12) months labor may be covered under this Limited Warranty. Refer to Owner's Manual for complete details.

**\*See Limited Warranty for complete details**

### ACCESSORIES:

**Deluxe Tech Chair:** Base: 24" diameter (61 cm) Height Range-Standard cylinder (excluding cushion): 14" to 16" (36 cm to 41 cm)  
Height Range: Higher cylinder-optional (excluding cushion): 19" to 24" (48 cm to 61 cm) Net weight: 15 lbs (7 kg).

**Shipping Dimensions:** 26" W (66 cm) x 26" L (66 cm) x 12" H (31 cm) Shipping Weight: 23 lbs (11 kg)

**Standard Tech Chair:** Base: 19" diameter (48 cm) Height Range-Standard cylinder (including cushion): 13" to 16" (33 cm to 41 cm)  
Height Range: Higher cylinder-optional (including cushion): 18" to 24" (46 cm to 61 cm)

Net weight: 14 lbs (6 kg).

**Shipping Dimensions:** 21" W (53 cm) x 21" L (53 cm) x 11" H (28 cm) Shipping Weight: 18 lbs (8 kg)

Note: Due to continuous improvements, all specifications are subject to change without notice. All dimensions and weights are approximate.